



# The first step into automation: Full software control of Stand-Alone DTs & DT Offline Systems

ERWEKA now offers a simple package to automate and future proof your existing stand-alone Dissolution Tester or Dissolution Offline System. The bundle comes with a Disso.NET 4.x preinstalled computer ready to connect to our DT or Offline System - with all important features activated: our powerful method and run editors, powerful user management and of course the 21 CFR part 11 conform Audit Trail and its powerful Audit Trail viewer.



Full Audit Trail according to 21 CFR Part 11



MS SQL Database



Support for USP methods 1, 2, 5 and 6



User Management

This is your first step into automation - let's take it together!

## Art. No. 27139 – Disso.NET 4.x for Stand-Alone DTs

- Upgrade your Stand-Alone DT with Audit Trail and full Software control!
- Complete upgrade package with Software Disso.NET 4.x with Computer, keyboard, mouse & screen
- Compatible with DTs 72x, 82x, 141x and 162x (check firmware compatibility)
- Includes all advanced Disso.NET 4.x features for stand-alone DTs, such as the expansive Audit Trail and the Audit Trail viewer
- Works with ASS-8 sampling station

## Art. No. 27154 – Disso.NET 4.x extension module for FRL

- Upgrade your existing DT Offline Systems with Disso.NET and Audit Trail and automate your Offline System!
- Works with FRL 82x and FRL 85x!
- requires Art. No. 27139

AltNumber	User	Time	UTC	Text	Parent	Child	Method	Details
2750	Operator	2016-10-19 10:39:04	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Operator	
2754	Operator	2016-10-19 10:39:04	+02:00:00	logged	Sequence 2016-09-24-HeadL	Run 2	Method	
2755	Operator	2016-10-19 10:40:24	+02:00:00	Status	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2756	Operator	2016-10-19 10:40:34	+02:00:00	Status	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2757	Operator	2016-10-19 10:40:34	+02:00:00	logged	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2758	Operator	2016-10-19 10:40:41	+02:00:00	logged	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2759	Operator	2016-10-19 10:40:48	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Operator	
2760	Administrator	2016-10-19 10:46:09	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Administrator	
2761	Administrator	2016-10-19 10:47:00	+02:00:00	Status	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2762	Administrator	2016-10-19 10:47:42	+02:00:00	changed	MeasuringSoftware	Dissolution Job	AllowTolerance	0.5
2763	Operator	2016-10-19 10:50:01	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Operator	
2764	Operator	2016-10-19 10:50:41	+02:00:00	Status	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2765	Operator	2016-10-19 10:50:50	+02:00:00	logged	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2766	Operator	2016-10-19 10:50:52	+02:00:00	Status	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2767	Operator	2016-10-19 10:50:52	+02:00:00	logged	Sequence 2016-09-24-HeadL	Run 2	Run	std646
2768	Operator	2016-10-19 10:50:52	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Operator	
2769	Administrator	2016-10-19 11:01:42	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Administrator	
2770	Administrator	2016-10-19 11:03:06	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Administrator	
2771	Operator	2016-10-19 11:03:59	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Operator	
2772	Operator	2016-10-19 11:05:26	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Operator	
2773	Administrator	2016-10-19 11:09:21	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Administrator	
2774	Administrator	2016-10-19 11:10:21	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Administrator	
2775	Operator	2016-10-19 11:13:27	+02:00:00	login	SettingMeasuringStation	SettingMeasuringStation	Operator	
2776	Operator	2016-10-19 11:13:27	+02:00:00	logged	SettingMeasuringStation	SettingMeasuringStation	Operator	

## Audit Trail Viewer

With the help of the Audit Trail Viewer, an authorized user can easily track all changes done to any audited data on the system.

## Dissolution Editor

The UI of our industry-proven Dissolution Editor has been refined and now offers a clearer focus for the task at hand.

